



MiOXSYS® System

Test	MiOXSYS®	
Catalog #	100283	
Technology	Electrochemical	
Target	Oxidation-Reduction Potential	
On Demand	Yes	
Testing Capacity	1	
Sample Type	Fresh	Frozen
	Liquefied Semen; RT	Liquefied Semen; Store at -80C
Sample Prep	Liquefied Semen; RT	Test Upon Thaw
TAT	5 minutes	
Hands on Time	< 1 minute	
Controls: Procedural	External Procedural Controls	
Sensitivity	Semen Parameters	
	63.4%	
Specificity	Semen Parameters	
	88.3%	
Sample Storage	Fresh	Frozen
	Liquefied Semen; RT	Liquefied Semen; Stored at -80C
Assay Storage	Room Temperature	



373 Inverness Parkway
Suite 206
Englewood, CO 80112

www.mioxsys.com

4350 Lassiter At North Hills
Suite 235
Raleigh, NC 27609



FOR EXPORT ONLY. NOT FOR SALE IN US.